



LTV 2 VENTILATOR SERIES

The ventilation Tools you Need throughout The Continuum of Care for Patients 10 kg and Above

With enhanced features, the new LTV2 ventilator series provides the power and portability needed for your early mobility and patient transport needs. The new LTV2 ventilators also provide improved patient synchrony in a compact, lightweight size with longer lasting hot swappable battery power.



Compact, Deployable Design

The small, lightweight LTV2 provides sophisticated critical care ventilation in Healthcare institutional use and intra-hospital transport environments. The LTV2 supports both invasive and non-invasive applications. It combines portability, impact-resistant construction, easy-to-read displays and comprehensive monitoring. A single solution for any patient weighing 10 kg or more, the LTV2 weighs only 11.5 lbs and is a pick-up-and-go solution that can travel with patients from field care to transport to the intensive care unit (ICU). It can take advantage of any 12-volt power source from cigarette lighters to outlets—and its internal battery allows it to be transferred from one power source to the next without interruption. Intuitive controls allow first responders to deploy the system with minimal retraining.



Repurpose ease

While the LTV2 can be stockpiled for use in emergencies, many of our customers use the LTV2 in hospital settings. More versatile than transport ventilators and significantly smaller than ICU systems, the LTV2 has intuitive controls and comprehensive monitoring capabilities that can play an essential role in day-to-day healthcare delivery and mass casualty medical care.

Oxygen (O2) resource management

The O2 Conserve feature reduces O2 consumption when the time and place of care require it. This feature helps extend the sometimes scarce resource of O2 during transport or emergencies and may be turned off when O2 resources are again in full supply.

Quick Patient Ventilation

Choose pediatric or adult to quickly configure initial ventilation before optimizing settings for each patient's needs. This feature aids in simple and effective emergency set-up.

PT MEDITRANS GLOBAL

Jakarta : Zamzam Square Unit R-5, Jl. Raya Condet No.4,
Jakarta 13530
Surabaya : Rungkut Menanggal Harapan Blok O No. 28,
Surabaya 60293
Other : Bandung, Makassar, Medan, Balikpapan

Hotline:
(021) 2287 9398
(031) 87856576
+62 812 8000 1102



Improved Portability and Enhanced Power Capabilities Make the LTV2 Ventilator Series the Right Choice for Any Situation

With enhanced features, the new LTV2 ventilator series provides the power and portability needed for your patient transport needs. The new LTV2 ventilators also provide improved patient synchrony in a compact, lightweight size with unlimited hot swappable battery power.

Hot Swappable Battery

- Uninterrupted freedom with 4-hour hot swappable battery
- 3.5-hour internal battery
- Two hot swappable batteries can simultaneously recharge in 2 hours
- Easy to read and accessible battery level check

LTV Series Accessories

Stand LTV Series Floor



Pediatric Circuit Without PEEP SPU



Heated Humidifier



Hepa Filter



Stand LTV Series Table



Adult Circuit Without PEEP SPU



PT MEDITRANS GLOBAL

Jakarta : Zamzam Square Unit R-5, Jl. Raya Condet No.4,
Jakarta 13530
Surabaya : Rungkut Menanggal Harapan Blok O No. 28,
Surabaya 60293
Other : Bandung, Makassar, Medan, Balikpapan

Hotline:
(021) 2287 9398
(031) 87856576
+62 812 8000 1102



KEMENKES RI AKL: 20403120816

PT. Meditrans Global

Jakarta : Zamzam Square Unit R-5, Jl. Raya Condet No.4,
Jakarta 13530

Surabaya : Rungkut Menanggal Harapan Blok O No. 28, Surabaya
60293

Other : Bandung, Makassar, Medan, Balikpapan

Hotline :

(021) 2287 9398

(031) 87856576

+62 812 8000 1102

Email :

info@meditransglobal.co.id