

# PROFEEL<sup>®</sup>

## LIGHTLY POWDERED

NR LATEX SURGICAL GLOVES

Powdered Latex Sterile Surgical Gloves

PROFEEL Lightly Powdered

COMFIT Powdered



Natural Rubber Latex

Lightly Powdered

Superior Sensitivity



Your Partner In Protection™

**Exceptional Dexterity**

Our innovative approach taken in meeting the healthcare professionals' increasingly stringent quality requirements is reflected in the outstanding features of **PROFEEL® Lightly Powdered NR Latex Surgical Gloves**.

Manufactured using specially designed porcelain formers, **PROFEEL® Lightly Powdered NR Latex Surgical Gloves** provide exceptional dexterity to enhance comfort and flexibility. Our proprietary manufacturing process ensures that the extractable NR latex protein content is kept to the minimum thereby reducing allergic skin reactions. In addition, the gloves contain a moderate level of dusting powder which facilitates swift removal of powder prior to surgical procedures.

• Textured surface on palm and fingers provide exceptional instrument grip



• Optimum thickness provides superior tactile sensitivity

• High tensile strength minimizes tear

• Reinforced cuff prevents tearing while donning

• Printed size and logo for easy identification

• Beaded cuff ensures easy donning



**PRODUCT INFORMATION\***

**Type:** Powdered and Gamma sterilized Surgical Gloves

**Primary Material:** Natural rubber latex

**Dimensions:**

Glove Size	Reorder #	Palm Width (mm)	Length (mm)
5½	P3055-20	73	281
6	P3060-20	78	282
6½	P3065-20	84	298
7	P3070-20	90	290
7½	P3075-20	96	298
8	P3080-20	104	300
8½	P3085-20	109	302
9	P3090-20	114	302

**Thickness:**

Location of Thickness Measurements	Single Wall (mm)
Finger (at 15 mm from the extreme tip)	0.23
Palm (at center of palm)	0.21
Cuff (at 25 mm from cuff end)	0.18

**Physical Properties:**

Parameters	Before Aging	After Aging
Tensile Strength (MPa)	28	22
Elongation (%)	860	820
Modulus 500% (MPa)	4	2.8
Force At Break (N)	18	16

**Powder:** US FDA approved USP absorbable cornstarch

**Color:** Natural (i.e. off-white)

**Design and Feature:** Hand specific, anatomically curved fingers, microroughed surface finish and beaded cuff (SH3)

**Packing:** 1 pair (1 left and 1 right) of gloves per inner wrap, 1 inner wrap per pouch, 50 pouches per dispenser, 6 dispensers per carton

**Pre-shipment Quality Inspection:**

Criteria	Insp. Level	AQL
Dimension	N=13 (EN455-2); S-2	4.0 (ASTM D3577)
Physical Properties	N=13 (EN455-2); S-2	4.0 (ASTM D3577)
1000ml Water Leak Test	G-1	1.0
Protein Content	N=8 (EN455-3)	N=3 (ASTM D3577)
Powder Content	N=2	N/A
Major Visual Inspection	G-1	1.5
Minor Visual Inspection	G-1	4.0

**Product Conformance:** In compliance with European Medical Device Directive 93/42/EEC (CE Class IIa) ASTM D3577

**Quality Assurance:** Manufacturing process is in compliance with US FDA Quality System Regulation (QSR) and BS EN ISO9001, BS EN ISO13485 Quality System

\*The product information provided is a guideline of typical performance values and characteristics of the product and not to be used as actual product specifications.



**WRP Asia Pacific Sdn Bhd** 1 4 7 8 1 7 V  
 Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi, 43900 Sepang, Selangor Darul Ehsan, MALAYSIA  
 Office +60-3-8706 1486  
 Facsimile +60-3-8706 1557  
 Email customer\_wrp@wrpworld.com



**REMESCO Handelges.m.b.H,**  
 Grinzinger Allee 5/25, 1190 Vienna, AUSTRIA  
 Office +43-1-328-50-88  
 Facsimile +43-1-328 29-49  
 Email office@remesco.com

Distributed by:

