

The Adjustable Aneroid Sphygmomanometer



ABN™

CALIBRA Plus

PRODUCT FEATURES

- Adjustable manometer that can easily reset the gauge back to zero needed
- Mini screwdriver included for stem of the gauge
- Complete with black nylon cuff adult and child size with heavy duty Latex inflation bladder and latex bulb
- Packed into black vinyl zipper case

PT SUGIH INSTRUMENTO ABADI
JL. TEMBOKAN RT. 1/01, DS. CIPEUNDEUY
PADALARANG 40553, INDONESIA

DISTRIBUTOR:
PT ABADINUSA USAHASEMESTA
JL. RADEN SALEH NO. 45G
JAKARTA 10330, INDONESIA

www.abadinusa.co.id
www.abnmedical.com