

PT. SUGIH	SPECIFICATIONS – ABN REGAL PLUS MOBILE CLOCK ANEROID SPHYGMOMANOMETER	
	Document No. : SIA-QA-13.00.399	Revision No. : 0
	Issue No. : 1	Revised Date : April 1, 2021
	Issued Date : April 1, 2021	Page : 1 of 2
	Arranged by : QA Staff	Signature :
	Approved by : QA Manager	Signature :

This specification is applicable for the following product types:

1. AS-081RP-M
2. AS-081RP-M-LF



- Clock gauge
- Bladder 2 tubes with nylon cuff
- Extended tubing with metal connector
- Cuff basket
- Coiled tube
- Bulb with deluxe valves
- Trolley

Component Specification			
No.	Description	Qty	Requirements
1	Clock aneroid gauge	1	<ul style="list-style-type: none"> <li>- Accuracy : <math>\pm 3</math> mmHg</li> <li>- Scale : <ul style="list-style-type: none"> <li>• 0 - 300 mmHg with 2 mmHg minor graduation</li> <li>• Mayor graduation at each 10 mmHg</li> <li>• Arabic numerals at each 20 mmHg division</li> </ul> </li> <li>- Durability : <math>\pm 3</math> mmHg throughout the pressure range after 10,000 alternating pressure cycles</li> <li>- Material : ABS</li> <li>- Color : Pantone cool gray 1C</li> <li>- Internal Mechanism : No pin stop</li> <li>- Faceplate : Black with luminescence scale</li> <li>- Lens material : Clear plastic</li> <li>- Serial number is imprinted on faceplate</li> </ul>
2	Cuff	2	<ul style="list-style-type: none"> <li>- Material: Nylon Oxford 210 D</li> <li>- Hook fastener tape material: Nylon</li> <li>- Loop fastener tape material: Nylon</li> <li>- Size: Adult and child</li> <li>- Arm circumference : 254 – 406 mm (adult), 184 – 267 mm (child)</li> <li>- Drawing: SIA-QA-13.00.192 (adult), SIA-QA-13.00.191 (child)</li> </ul>

<b>PT. SUGIH</b>	<b>SPECIFICATIONS – ABN REGAL PLUS MOBILE CLOCK ANEROID SPHYGMOMANOMETER</b>	
	Document No. : SIA-QA-13.00.399	Revision No. : 0
	Issue No. : 1	Revised Date : April 1, 2021
	Issued Date : April 1, 2021	Page : 2 of 2
	Arranged by : QA Staff	Signature :
	Approved by : QA Manager	Signature :

3	Bladder, 2 tubes	2	- Material : Natural latex <u>OR</u> polychloroprene latex - Color : Black <u>OR</u> terracotta - Drawing: SIA-QA-13.00.004 (Latex adult), SIA-QA-13.00.003 (Latex child), SIA-QA-13.00.009 (LF adult), SIA-QA-13.00.008 (LF child)
4	Inflation bulb	1	- Material : Natural latex <u>OR</u> polychloroprene latex - Color : Black <u>OR</u> terracotta - Drawing: SIA-QA-13.00.268 (Latex), SIA-QA-13.00.018 (LF)
5	Coiled tube	1	- Material : Natural latex <u>OR</u> TPU - Color : Black - Drawing: SIA-QA-13.00.036 (Latex large 8ft), SIA-QA-13.00.173 (TPU)
6	End valve	1	- Material : Ni chrome-plated brass - Drawing : SIA-QA-13.00.020 - Type : Deluxe
7	Air release valve	1	- Material : Ni chrome-plated brass - Drawing : SIA-QA-13.00.022 - Type : Deluxe
8	Connector	2 sets	- Material : Brass rod - Finishing : Ni chrome plating - Drawing : SIA-QA-13.00.303
9	Extended tube	1	- Material : Natural latex <u>OR</u> polychloroprene latex - Color : Black - Length : 10 cm
10	Cuff basket	1	- Material : ABS - Color : Pantone cool gray 1C
11	Trolley	1 set	- Material : STKM, stainless steel - Color : Pantone cool gray 1C
12	Instruction Manual	1	- Material : Artpaper - Drawing: SIA-QA-13.00.281
13	Warranty Card	1	- Material : BC - Drawing: SIA-QA-13.00.297
14	Carton	1	- Material : 2-ply corrugated paper - Drawing: SIA-QA-13.00.392