

ABN™

REGAL Plus

CLOCK ANEROID SPHYGMOMANOMETER

The Stylish Design Large Clock
Aneorid Sphymomanometer
for Professional



- Featuring large face 8", 300 mmHg no pin stop mechanism of manometer gauge
- The large dial plate for easy reading supported with the attractive strong base design, make this model grows to be the favourite for hospitals or clinics
- Comes with complete inflation system of the superior quality Adult and Child black nylon cuff with heavy duty latex inflation bladder, bulb and coiled tubing, they can all be placed nicely into storage compartment
- Designed as the desk model to get the flexibility of placement in all operations, Wall model to be attached to 90° swivel wall mounting bracket and Mobile model The Moveable Clock Aneroid Sphymomanometer

PT SUGIH INSTRUMENTO ABADI
JL. TEMBOKAN RT. 1/01, DS. CIPEUNDEUY
PADALARANG 40553, INDONESIA

DISTRIBUTOR:
PT ABADINUSA USAHASEMESTA
JL. RADEN SALEH NO. 45G
JAKARTA 10330, INDONESIA

www.abadinusa.co.id
www.abnmedical.com