

PT. SUGIH	SPECIFICATIONS – ABN CALIBRA PLUS ANEROID SPHYGMOMANOMETER	
	Document No. : SIA-QA-13.00.395	Revision No. : 0
	Issue No. : 1	Revised Date : April 1, 2021
	Issued Date : April 1, 2021	Page : 1 of 2
	Arranged by : QA Staff	Signature :
	Approved by : QA Manager	Signature :



- Latex bladder 2-tube with black nylon cuff
- Gauge
- Latex bulb with standard valves
- Mini screwdriver

Component Specification			
No.	Description	Qty	Requirements
1	Gauge	1	<ul style="list-style-type: none"> - Accuracy: ± 3 mmHg - Scale: <ul style="list-style-type: none"> • 0-300 mmHg with 2 mmHg minor graduation. • Major graduation at each 10 mmHg. • Arabic numerals at each 20 mmHg division. - Durability: ± 3 mmHg throughout the pressure range after 10.000 alternating pressure cycles. - Material: Zinc Alloy - Coating/Finish: Powder coating, black. - Internal Mechanism: No stop Pin. - Lens Material: Clear plastic - Gauge clip: Stainless steel
2	Cuff	2	<ul style="list-style-type: none"> - Material: Nylon Oxford 210 D - Hook fastener tape material: Nylon - Loop fastener tape material: Nylon - Size: Adult and child - Arm circumference : 254 – 406 mm (adult), 184 – 267 mm (child) - Drawing: SIA-QA-13.00.192 (adult), SIA-QA-13.00.191 (child)
3	Latex Bladder 2-tube	2	<ul style="list-style-type: none"> - Material: Natural Latex - Color: Black - Size: Adult and child - Drawing: SIA-QA-13.00.004 (adult), SIA-QA-13.00.003 (child)
4	Inflation Bulb	1	<ul style="list-style-type: none"> - Material : Natural Latex - Warna : Hitam - Drawing : SIA-QA-13.00.268
5	End Valve	1	<ul style="list-style-type: none"> - Material: Ni chrome-plated brass - Type: Standard - Drawing: SIA-QA-13.00.019
6	Air Release Valve	1	<ul style="list-style-type: none"> - Material: Ni chrome-plated brass - Type: Standard - Drawing: SIA-QA-13.00.021
7	Mini screwdriver	1	<ul style="list-style-type: none"> - Material: Steel & plastic

PT. SUGIH		SPECIFICATIONS – ABN CALIBRA PLUS ANEROID SPHYGMOMANOMETER	
		Document No. : SIA-QA-13.00.395	Revision No. : 0
		Issue No. : 1	Revised Date : April 1, 2021
		Issued Date : April 1, 2021	Page : 2 of 2
		Arranged by : QA Staff	Signature :
		Approved by : QA Manager	Signature :

8	Zippercase	1	- Material: Vinyl leather - Color: Black - Drawing: SIA-QA-13.00.012
9	Warranty Card	1	- Material: BC - Drawing: SIA-QA-13.00.297
10	Instruction Manual	1	- Material: Artpaper - Drawing: SIA-QA-13.00.064
11	Packaging	1	- Material: E-flute corrugated with duplex paper - Drawing: SIA-QA-13.00.390