

PT. SUGIH	SPECIFICATIONS – ABN DM-500 PLUS DESK DIGITAL SPHYGMOMANOMETER	
	Document No. : SIA-QA-13.00.396	Revision No. : 0
	Issue No. : 1	Revised Date : April 1, 2021
	Issued Date : April 1, 2021	Page : 1 of 2
	Arranged by : QA Staff	Signature :
	Approved by : QA Manager	Signature :



Digital manometer

Cuff & bladder

Bulb with valve assy

Cuff compartment

Specifications			
No.	Component	Qty.	Description
1	Digital Manometer with Cuff Compartment	1 set	<ul style="list-style-type: none"> - Material : ABS Model : DM-500P Indicator Column : LCD Screen Scale : 20 – 280 mmHg Accuracy : ± 3 mmHg Inflation : Bulb w/ valve Deflation : Air Release Valve Power Supply : AA-size Batteries (2 pcs) Power Consumption : 0.5 W (max.) Voltage : DC 3V Operating Environments : +10°C s/d +40°C, 85% humidity Storage Environments : -5°C s/d +50°C, 85% humidity Weight : 780 grams (excl. battery)
2	PVC Tubing	1	<ul style="list-style-type: none"> - Material : PVC - Color : Black - Length : 25 cm - ID : 4 mm - OD : 7.5 mm
3	Air Release Valve	1	<ul style="list-style-type: none"> - Material : Ni chrome-plated brass - Drawing : SIA-QA-13.00.022 - Type : Deluxe - With ABN Stamp
4	End Valve	1	<ul style="list-style-type: none"> - Material : Ni chrome plated brass - Drawing : SIA-QA-13.00.020 - Type : Deluxe

PT. SUGIH		SPECIFICATIONS – ABN DM-500 PLUS DESK DIGITAL SPHYGMOMANOMETER	
		Document No. : SIA-QA-13.00.396	Revision No. : 0
		Issue No. : 1	Revised Date : April 1, 2021
		Issued Date : April 1, 2021	Page : 2 of 2
		Arranged by : QA Staff	Signature :
		Approved by : QA Manager	Signature :

5	Cuff	2	<ul style="list-style-type: none"> - Material: Nylon Oxford 210 D - Hook fastener tape material: Nylon - Loop fastener tape material: Nylon - Size: Adult and child - Arm circumference : 254 – 406 mm (adult), 184 – 267 mm (child) - Drawing: SIA-QA-13.00.192 (adult), SIA-QA-13.00.191 (child)
6	PVC Bladder 2-tube	2	<ul style="list-style-type: none"> - Material : PVC (with black PVC Tubing 50 cm) - Size : Adult and child - Reference : SIA-QA-13.00.243 (adult), SIA-QA-13.00.242 (child)
7	Bulb	1	<ul style="list-style-type: none"> - Material : Natural Latex - Warna : Hitam - Drawing : SIA-QA-13.00.268
8	Connector	1 set	<ul style="list-style-type: none"> - Material : Brass Rod - Finishing : Ni-chrome plating - Drawing: SIA-QA-13.00.303
9	Battery	2 pcs	<ul style="list-style-type: none"> - Type : AA Battery
10	Instruction Manual	1	<ul style="list-style-type: none"> - Material : HVS - Bilingual version (Bahasa & English) - Drawing: SIA-QA-13.00.196
11	Warranty Card	1	<ul style="list-style-type: none"> - Material : BC - Drawing: SIA-QA-13.00.297
12	Unit Box	1	<ul style="list-style-type: none"> - Material : E-flute with duplex cover - Drawing: SIA-QA-13.00.393